

NASAL CANNULAS

OXYGEN-CONSERVING CANNULAS

Description:

The Hudson RCI® Oxygen-Conserving Cannula incorporates a patented split-prong nasal design and two oxygen supply tubes with standard connectors.

- Design provides control and comfort in synchronizing breathing with oxygen conserving regulators
- Use of this cannula requires a conserving regulator with two ports, one for oxygen delivery and one for pressure sampling
- Individually packaged



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Tubing Length	Standard Tubing	Qty/Case
1850	5 ft (2)	•	25
1851	7 ft (2)	•	25

© 2021 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.