

SMALL VOLUME NEBULIZERS

Teleflex offers a suite of Hudson RCI® Small Volume Nebulizers that promote positive patient outcomes and help to maximize clinician efficiency.



ISO-NEB® FILTERED NEBULIZER SYSTEM

The ISO-Neb® Filtered Nebulizer System is designed to reduce environmental contamination during aerosol therapy with inspiratory and expiratory flows that are isolated by one-way valves with exhalation directed through an efficient filter system.

[View Details](#)

MICRO MIST® NEBULIZER

The Micro Mist® Nebulizer is a high-performing nebulizer that delivers an optimal particle size at angles up to 90°. The unique design maximizes options for placement, without sacrificing aerosol quality. It is used for intermittent aerosol therapy and can be relied upon to deliver both performance and economy.

[View Details](#)

UP-DRAFT II® OPTI-NEB® NEBULIZER

Designed to be a universal nebulizer, the Up-Draft II® Opti-Neb® Nebulizer gives you added flexibility in flow rates, delivering quality aerosols at flows ranging from 5 to 9 LPM and making it a popular choice for many facilities looking to consolidate options into one high-performance product.

[View Details](#)

UP-DRAFT® NEBULIZER

The Up-Draft® Nebulizer is a high-output nebulizer available in a variety of configurations.

[View Details](#)