

NASAL CANNULAS

SOFTECH® PLUS ETCO₂ SAMPLING CANNULA

Description:

The Softech® Plus ETCO₂ Sampling Cannula allows delivery of supplemental oxygen while simultaneously sampling expired CO₂, and features softer nasal prongs and lariat for enhanced patient comfort. The split nare design allows separation of the CO₂ readings and the delivery of oxygen, and helps produce a sharp wave form for clinician diagnostics.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat Am

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
2842	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 2 in. CO ₂	Male	50
2843	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 2 in. CO ₂	Female	50
2844	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Female	50
2845	Softech® Plus ETCO ₂ Sampling Cannula, Adult	7 ft O ₂ / 7 ft CO ₂	Male	50
2846	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 2 in. CO ₂	Female	50
2847	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 14 ft CO ₂	Female	50

Catalog #	Size	Tubing Length	Luer Connector	Qty/Case
2848	Softech® Plus ETCO ₂ Sampling Cannula, Adult	14 ft O ₂ / 14 ft CO ₂	Male	50
2849	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 2 in. CO ₂	Female	50
2850	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 7 ft CO ₂	Female	50
2851	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	7 ft O ₂ / 7 ft CO ₂	Male	50
2852	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 2 in. CO ₂	Female	50
2853	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 14 ft CO ₂	Female	50
2854	Softech® Plus ETCO ₂ Sampling Cannula, Pediatric	14 ft O ₂ / 14 ft CO ₂	Male	50
2855	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 2 in. CO ₂	Female	50
2856	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 7 ft CO ₂	Female	50
2857	Softech® Plus ETCO ₂ Sampling Cannula, Infant	7 ft O ₂ / 7 ft CO ₂	Male	50

© 2019 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions; contact customer service to confirm availability in your region.