

CATALOG # : 2850

## RESPIRATORY

### Description:

Softech® Plus ETCO<sub>2</sub> Sampling Cannula, Pediatric, 50/CS (Individually Packaged)



United States

To order this product, please call:

**866-246-6990**

Lat

Am

### Details:

- Pediatric Cannula with 7 ft oxygen supply tubing and 7 ft CO<sub>2</sub> sampling line with a female luer connector

© 2021 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.