ENDOTRACHEAL TUBE ACCESSORIES

RUSCH® SILKOSPRAY™ UNIVERSAL SILICONE SPRAY

Description:

The Rusch® Silkospray™ Universal Silicone Spray reduces friction on the surface of devices made of rubber, latex and PVC. Silkospray is CFC-free, and is not intended for use with devices made of silicone.



To order this product, please call:

+1 (919) 433-4999

Catalog #	Description	Qty/Box
556000-000000	Rusch® Silkospray™ Universal Silicone Spray, 500 mL	10

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.